

ASTM D 482 - D 874

SCOPE:

This test method covers the determination of ash in the range 0.001-0.180 mass %, from distillate and residual fuels, gas turbine fuels, crudes oils, lubricating oils, waxes, and other petroleum products, in which any ash-forming materials present are normally considered to be undesirable impurities or contaminants.

**ASH FROM PETROLEUM PRODUCTS
- REF 23963 -**

MAIN CHARACTERISTICS

Muffle furnace model L3/11, 3 liters capacity, maximum temperature of 1100°C.

NECESSARY ACCESSORIES (for ASTM D482)

- REF 10532** Porcelain dish, 100 mL, flat bottom (pack of 6)
- REF 10687** Stainless curved crucible tong (pack of 2)
- REF 10702** Glass balls diam. 2mm (1kg)
- REF 10770** Porcelain plate for below dessicator diam. 200mm
- REF 9411300** Platinum crucible (100mL)

NECESSARY ACCESSORIES (for D874)

- REF 10508** Porcelain crucible, 40mL tall form (pack of 6)
- REF 10509** Porcelain crucible, 80mL tall form (pack of 6)
- REF 10510** Porcelain crucible, 120mL tall form (pack of 6)
- REF 10515** Porcelain crucible, 80mL low form (pack of 6)
- REF 10532** Porcelain dish, 100mL, flat bottom (pack of 6)
- REF 10687** Stainless curved crucible tong (pack of 2)
- REF 10702** Glass balls diam. 2mm (1kg)
- REF 10770** Porcelain plate for below dessicator
- REF 27203** Soda lime glass dessicator diam. 200mm
- REF 15377** Silica crucible (50 mL tall form)
- REF 15386** Silica crucible (75 mL tall form)
- REF 15387** Silica crucible (130 mL tall form)
- REF 15671** Platinum crucible (50mL)

SPARE PARTS

- REF 23964** Set of heating elements
- REF 23965** Thermocouple K (160mm)



NORMALAB FRANCE SAS
14 rue des Lilas - F-76210 Lintot
Tel. : +332 35 38 59 59 - Fax : +332 35 38 78 55

NORMALAB BELGIUM S.A.
Parc Industriel
28 rue des Pieds-d'Alouette – BE-5100 NANINNE / Belgium
Tel : +32 81 71 19 90 - Fax : +32 81 22 27 01

DISTRIBUTED BY

NORMALAB FRANCE SAS
14 rue des Lilas - F-76210 Lintot
Tel. : +332 35 38 59 59 - Fax : +332 35 38 78 55

NORMALAB BELGIUM S.A.
Parc Industriel
28 rue des Pieds-d'Alouette – BE-5100 NANINNE / Belgium
Tel : +32 81 71 19 90 - Fax : +32 81 22 27 01

DISTRIBUTED BY