



- ⊕ **Fast cleaning (5 .. 15 min)**
- ⊕ **6 positions**
- ⊕ **BIO degradable solvent**
- ⊕ **Negligible solvent consumption**

General

Designed in order to eliminate hand cleaning of glass kinematic viscometer tubes

Viscometer washer supplied with 6 nozzle stoppers. It allows external and internal washing of all current types of viscometer tubes. Viscometer tubes are suspended in BIOSANE T215 solvent vapors at 80°C

Ordering code

18450

Recommended solvent

25699 Biosane T 215 solvent (per can of 5 l)

Specs

For use on AC 230V, 50/60 Hz - 5A (115V on request)

Dimensions (wxdxh) 240x260x350 mm

Weight : ± 7 kg